



# Fully Automated Magtration System GC series for total RNA Purification from Animal Tissue

## Performance of Total RNA purification from mouse liver and kidney

Table 1 Materials and Reagents

|                     |  |
|---------------------|--|
| Reagent             | MagDEA RNA 100 Cell/Tissue (GC)<br>(Code No. E7004)  |
| Protocol            | MagDEA RNA 100 Cell/Tissue ver. 1.0  |
| Sample              | Mouse Liver and Kidney   |
| Sample Pretreatment | Homogenize with <i>RNAlater</i> <sup>®</sup>   |
| Sample Volume       | 100 $\mu$ l (Liver: 8.1 mg, Kidney: 7.3 mg)  |
| Elution Volume      | 50 $\mu$ l / 100 $\mu$ l   |
| Operation Time      | Approx. 40 min. (Without DNase I treatment)<br>Approx. 65 min. (Including DNase I treatment) |




Figure 1. Agarose gel electrophoresis of purified total RNA treated with DNase I from mouse liver and kidney by Magtration System 12GC

## Yields and Quality of purified total RNA

Table 2 Yield and Purity (N=3)

| Sample       | Yield ( $\mu$ g) | $A_{260}/A_{280}$ | $A_{260}/A_{230}$ |
|--------------|------------------|-------------------|-------------------|
| Mouse Liver  | 25.8 $\pm$ 0.5   | 2.05 $\pm$ 0.01   | 1.78 $\pm$ 0.03   |
| Mouse Kidney | 10.9 $\pm$ 0.3   | 2.06 $\pm$ 0.02   | 1.91 $\pm$ 0.03   |

Table 3 Quality of purified RNA ( Agilent BioAnalyzer 2100)

| Sample       | rRNA ratio [28s/18s] | RNA Integrity number |
|--------------|----------------------|----------------------|
| Mouse Liver  | 1.30 $\pm$ 0.00      | 8.5 $\pm$ 0.1        |
| Mouse Kidney | 1.26 $\pm$ 0.02      | 8.3 $\pm$ 0.2        |




Figure 2. Electropherogram of total RNA isolated from mouse liver and kidney by Agilent 2100 bioAnalyzer.

## RT-PCR amplification using total RNA purified from the tissue by Magtration System 12GC

Table 4 Condition of PCR amplification

| Template      | 2 $\mu$ l   | RT-PCR Cycle |
|---------------|---|--------------|
| Primer Target | RT- PCR of Mouse G3PDH gene (452 bp) using OneStep RT-PCR Kit (Qiagen Inc.) | RT step      |
|               |   | Repeat       |
|               |   | Extension    |
|               |   | 30 cycles    |
|               |   |              |




Figure3. RT-PCR of total RNAs purified from mouse liver and kidney .

### Advantages

The Magtration System GC series is the most reliable genomic DNA purification system in the current market. The Magtration System GC series enables processing of sample in a short time in a small footprint on a lab-bench. The genomic DNA purified by the Magtration System GC series is sufficient in yield and purity to be used in downstream applications directly, such as PCR. The Magtration System GC series is a true walk-away automation system.



Asia Region  
Precision System Science Co., Ltd.  
Tech-support@pss.co.jp  
URL: <http://www.pss.co.jp>



United States  
Precision System Science USA, Inc.  
contact@pssbio.com  
URL: <http://www.pssbio.com>



European Region  
Precision System Science Europe GmbH  
contact@pss-europe.de  
URL: <http://www.pss-europe.de>