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URL : <http://www.pss.co.jp>, E-mail : service@pss.co.jp

ENGLISH

This product is the Protocol IC Card for the automated nucleic acid extraction system magLEAD 6gC or magLEAD 12gC. Please surely read and understand this document and operation manual of the system prior to use.

Symbols;

: Caution

: Contents

: Manufacturer

: Refer to Instructions for use

: CE mark for Europe

: In Vitro Diagnostic medical devices (98/79/EC)

: Authorized Representative in the European Community

1. Abstract of product

This Protocol IC Card is for the nucleic acid extraction by the use of PSS produced pre-filled nucleic acid extraction reagent and the automated nucleic acid extraction system magLEAD 6gC or magLEAD 12gC. (As for the nucleic acid extraction reagent, please refer to the instructions for use of extraction reagent.)

2 How to use this product

2.1 Safety instructions

- This Protocol IC Card is dedicated for the automated nucleic acid extraction system magLEAD 6gC or magLEAD 12gC. Therefore, please read carefully the operation manual of the system prior to use.
- The Protocol IC Card is a dedicated for the instrument model each. It cannot be used to non-compliant model. Please make sure the instrument model described in the label prior to use.
- In case of instrument error messages, please refer to the instrument operation manual.
- This product contains electrical board. Therefore, please do not drop (giving a shock), bend, wet and expose to extreme temperature and sunlight.
- Please make sure the power switch of the instrument is OFF when inserting or removing the IC card.

2.2 Operation procedures

1. Insert Protocol IC Card to the instrument.

NOTE: Please make sure the power switch of the instrument is OFF.

2. Turn ON the instrument.
3. Please turn off the power switch of the instrument when removing the Protocol IC Card from the instrument

These explanations are based on the status of Feb. 2016. Please be noted information such as specification can be changed without information.

Manufacturer

Precision System Science Co., Ltd.

Main Office

Address Zip code: 271-0064

88, Kamihongou, Matsudo City, Chiba, Japan

Contact

TEL : +81-(0)47-303-4800, FAX : +81-(0)47-303-4811

URL : <http://www.pss.co.jp>, E-mail : service@pss.co.jp

Emergo Europe

Address

Molenstraat 15, 2513 BH The Hague, The Netherlands