Instructions for Use

magLEAD Consumable Kit

Version 1.1
December 2015

CONT

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ENGLISH

This product is the plastic consumable kit for the automated nucleic acid extraction system magLEAD 6gC and magLEAD 12gC PSS produced. Please surely read and understand this document and operation manual of the system prior to use.

Contents

1. Symbols	5
2. Abstract of product	5
2.1 Intended for Use	5
2.2. Components of product	5
2.3. Storage instructions	5
3 How to use this product	5
3.1. Safety instructions	5
3.2. Operation procedures	6

1. Symbols

: Caution

: Lot number

: Product code

: Contents

: Refer to Instructions
for use

: Do not reuse

: Expiration period

: Keep away from
sunlight

: Manufacturer

2. Abstract of product

2.1 Intended for Use

This plastic consumable kits is for the nucleic acid extraction by the use of PSS produced pre-filled nucleic acid extraction reagent "MagDEA® Dx SV" and automated extraction system magLEAD 6gC and magLEAD 12gC. (As for the nucleic acid extraction reagent, please refer to the instructions for use of extraction reagent.)

2.2. Components of product

- | | |
|---|---------------------------|
| 1. Tip and Sheath set | 52 sets |
| 2. Micro tube 1.5mL, PP (Electron irradiated)
For Sample tube and elution tube | 100 pcs (50 pcs × 2 bags) |
| 3. Screw Cap (Electron irradiated)
Screw cap for tubes above | 100 pcs (50 pcs × 2 bags) |

2.3. Storage instructions

Please store this product indoors with avoiding high-temperature and humidity. Additionally, please keep away from sunlight.

3 How to use this product

3.1. Safety instructions

Please surely read and understand the followings prior to usage.

- This plastic consumables are dedicated for the fully automated nucleic acid extraction system PSS produced. Therefore, please read carefully the operation manual of the system prior to use.
- In case of instrument error messages, please refer to the instrument operation manual.
- Please do not reuse the used consumables.
- Please do not use the deformed and/or discolored consumables. Furthermore, if you find an attached foreign materials, please do not use. If such consumables is used, there is a risk to affect the performance and operation/run.

NOTE

- Please discard the used consumables in accordance with a regulation of your lab and/or regions.

3.2. Operation procedures

Please read carefully the operation manual of the system prior to use.

1. Insert Protocol IC Card to the instrument, then, turn on the instrument.
2. Set the consumables (This product), reagent cartridges (MagDEA[®] Dx SV) and samples to the instrument, according to the message shown in the display.

E : Micro tube (for elution)
T1 : Tip and Sheath set
T2 : -
S : Micro tube (for sample)

3. The quantity needed for the one sample preparation is below.

Tip and Sheath set	1 set
Micro tube (for elution tube)	1 pc
Micro tube (for sample tube)	1 pc
MagDEA [®] Dx SV (reagent cartridge)	1 pc
4. Please make sure the all plastic consumables, reagent cartridges and samples are set correctly.
5. Close the front door of the instrument.
6. Press START button. The extraction run will be started.
7. After completion of the run, remove the elution tube, used MagDEA[®] Dx SV and other consumables.

Magtration[®] and MagDEA[®] are registries by Precision System Science Co., Ltd. These explanations are based on the status of Dec. 2015. Please be noted information such as specification can be changed without information.

Produced by / sold by

Precision System Science Co., Ltd.
88 Kamihongou Matsudo-shi, Chiba 271-0064 JAPAN
Tel: +81(0)47-303-4800 Fax: +81(0)47-303-4811
URL : <http://www.pss.co.jp>
E-mail : service@pss.co.jp

Precision System Science USA, Inc.
5673 West Las Positas Blvd., Suite 202, Pleasanton, CA 94588, U.S.A.
E-mail : contact@pssbio.com

Precision System Science Europe GmbH
Steinstrasse 27, D-40210 Düsseldorf, Germany
E-mail : contact@pss-europe.de